

# CryoSolve

## cGMP organic cryopreservation media


### RAISING THE STANDARDS OF A PROVEN FORMULA

#### WHY USE CRYOSOLVE?

- Single-Use syringes:
  - Eliminate sharps
  - Avoid cross contamination and maintain sterility
  - Save time
  - Streamline process steps
- All organic raw materials - serum-free/animal free



#### VIABILITY Average %



- Manufactured under cGMP – ensuring that you are using only the highest quality material
- Raw materials USP/EP grade
- All equipment passed compatibility studies
- Ready-to-Use\*
- Average 97% cell viability
- Made in USA

### The Latest Generation in Cell Preservation

#### SPECIFICATIONS

- CryoSolve Syringes filled with: DMSO USP/EP 55% w/v, Dextran-40 USP/EP 5% w/v and WFI
- 10% final DMSO concentration
- Ready-to-use syringes containing 7.0 mL CryoSolve
- Storage Temperature: 2°C to 8°C

\*Refer to Method of Use for instructions

